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1. Rosenquist C., Qvist P., Bjarnason N.H., and Christiansen C. 1995. Measurement of a more stable region of osteocalcin in serum by ELISA with two monoclonal antibodies. *Clin Chem* 41, 10:1439-1445
2. Takahashi M, Kushida K, Nagano A, Inoue T. 2000. Comparison of the analytical and clinical performance characteristics of an N-MID versus an intact osteocalcin immunoradiometric assay. *Clin Chim Acta*. 294:67-76
3. Marcus R; Holloway L; Wells B; Greendale G; James MK; Wasilaukas C; Kelaghan J. 1999. The Relationship of Biochemical Markers of Bone Turnover to Bone Density Changes in postmenopausal Women: Results from the Postmenopausal Estrogen/Progestin Interventions (PEPI) Trial. *Journal of Bone and Mineral Research*; 14 p1583-1595
4. McClung M; Clemmesen B; Daifotis A; Gilchrist NL; Eisman JA; Weinstein RS; Fuleihan Ge; Reda C; Yates AJ; Ravn P. 1998. Alendronate Prevents Postmenopausal Bone Loss in Women Without Osteoporosis. A Double-Blind, Randomized, Controlled Trial. Alendronate Osteoporosis Prevention Study Group. *Annals of Internal Medicine*; 128(4) p253-261
5. Meunier P; Confravreux E; Tupinon I; Hardouin C; Delmas PD; Balena R. 1997. Prevention of Early Postmenopausal Bone Loss with Cyclical Etidronate Therapy (A Double-Blind, Placebo-Controlled Study and 1-Year Follow-Up). *Journal of Clinical Endocrinology & Metabolism*; 82 p2784-2791
6. N. H. Bjarnason, K. Bjarnason, J. Haarbo, C. Rosenquist, and C. Christiansen. 1996. Tibolone: Prevention of Bone Loss in Late Postmenopausal Women. *J Clin Endocrinol Metab* 81:2419-2422
7. Ravn P; Clemmesen B; Riis BJ; Christiansen C. 1996. The Effect on Bone Mass and Bone Markers of Different Doses of Ibandronate: A New

8. Bisphosphonate for Prevention and Treatment of Postmenopausal Osteoporosis: A 1-year, Randomized, Double-blind, Placebo-controlled Dose-finding Study. *Bone* 19 p527-533
8. Christiansen, C., L.B. Tanko, L. Warming, A. Moelgaard, S. Christgau, P. Qvist, M. Baumann, L. Wiczorek, and N. Hoyle. 2003. Dose dependent effects on bone resorption and formation of intermittently administered intravenous ibandronate. *Osteoporos Int*
9. Junqueira PA, da Fonseca AM, Bagnoli VR, Giannella-Netot D, Manguiera CL, Coimbra CN, Pinotti JA. 2002. Comparison of bone remodeling indicators in climacteric women. *Int J Fertil Womens Med* 47:174-181
10. Ohishi T, Takahashi M, Kushida K, Hoshino H, Tsuchikawa T, Naitoh K, Inoue T. 1998. Changes of biochemical markers during fracture healing. *Arch Orthop Trauma Surg*. 118:126-130
11. Ravn P, Christensen JO, Baumann M, Clemmesen B. 1998. Changes in biochemical markers and bone mass after withdrawal of ibandronate treatment: prediction of bone mass changes during treatment. *Bone*. 22:559-564
12. Ravn P, Clemmesen B, Christiansen C. 1999. Biochemical markers can predict the response in bone mass during alendronate treatment in early postmenopausal women. Alendronate Osteoporosis Prevention Study Group. *Bone* 24:237-244
13. Ravn P; Rix M; Andreassen H; Clemmesen B; Bidstrup M; Gunnes M. 1997. High bone turnover is associated with low bone mass and spinal fracture in postmenopausal women. *Calcif Tissue Int*. 60:255-260
14. Tanko, L.B., U. Mouritzen, H.J. Lehmann, L. Warming, A. Moelgaard, S. Christgau, P. Qvist, M. Baumann, L. Wiczorek, N. Hoyle, and C. Christiansen. 2003. Oral ibandronate: changes in markers of bone turnover during adequately dosed continuous and weekly therapy and during different suboptimally dosed treatment regimens. *Bone* 32:687-693

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生产者名称: Immunodiagnostic Systems Limited (IDS Ltd)
生产者/生产场所地址: 10 Didcot Way, Boldon Business Park, Boldon, Tyne & Wear, NE35 9PD, UK

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