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1. Ceunick F De, Sabatini M, Renoux V, Nanteuil G de, Pastoureau P. Urinary collagen type II C-telopeptide fragments are sensitive markers of matrix metallo-proteinase dependent cartilage degradation in rat adjuvant induced arthritis. *J Rheumatol* (2003); 30: 1561-1564.
2. Christgau S, Tankó LB, Cloos PAC, Mouritzen U, Christiansen C, Delaissé J-M, Høegh-Andersen P. Suppression of Elevated Cartilage Turnover in Postmenopausal Women and in Ovariectomized Rats by Estrogen and a Selective Estrogen Receptor Modulator (SERM). Submitted.
3. Christgau S, Henrotin Y, Henriksen DB, Rovati LC, Collette J, Bruyere O, Deroisy R, Christiansen C, Reginster JY. Cartilage Degradation In Glucosamine Sulphate Treated Knee Osteoarthritis Patients With Elevated Levels Of Urinary Collagen Type II Ctelopeptide Fragments. *Clin Exp Rheumatol*. 2004 Jan-Feb;22(1):36-42
4. Christgau S, Garnero P, Fledelius C, Moniz C, Rosenquist C, Ensig M, Gineyts E, Christiansen C, Qvist P. Collagen type II degradation products in urine as an index of cartilage degradation. *Bone* (2001); 29: 209-215.
5. Forsblad d'Elia H, Christgau S, Mattsson L-Å, Saxne T, Ohlsson C, Nordborg E, Carlsten H. Hormone replacement therapy decreases markers of cartilage and bone metabolism in rheumatoid arthritis. Submitted.
6. Garnero P, Mazières B, Gueguen A, Abbal M, Berdah L, Freiburghaus C, Lequesne M, Nguyen M, Salles J-P, Vignon E, Dougados M. Association of 10 molecular markers of bone, cartilage and synovium with disease activity and joint damage in hip osteoarthritis patients: the ECHODIAH cohort. *ACR* 2003.
7. Garnero P, Landewé R, Boers M, Verhoeven A, Linden S van der, Christgau S, Heijde D van der, Boonen A, Geusens P. Association of baseline levels of markers of bone and cartilage degradation with long-term progression of joint damage in patients with early rheumatoid arthritis: the COBRA Study. *Arthritis & Rheum* (2002); 46:2847-2856.
8. Garnero P, Ayrat X, Rousseau J-C, Christgau S, Sandell L, Delmas PD, Dougados M. Uncoupling of type II collagen synthesis and degradation predicts progression of joint damage in patients with knee osteoarthritis. *Arthritis & Rheum* (2002); 46:2613-2624.
9. Garnero P, Gineyts E, Christgau S, Finck B, Delmas PD. Association of baseline levels of urinary glucosyl-galactosyl pyridinoline and type II collagen C-telopeptide with progression of joint destruction in patients with early rheumatoid arthritis. *Arthritis & Rheum* (2002); 46: 21-30.
10. Garnero P, Christgau S, Delmas PD. The bisphosphonate Zoledronate decreases type II collagen breakdown in patients with Paget's disease of bone. *Bone* (2001); 28: 461-464.
11. Garnero P, Piperno M, Gineyts E, Christgau S, Delmas PD, Vignon E. Cross sectional evaluation of biochemical markers of bone, cartilage, and synovial tissue metabolism in patients with knee osteoarthritis: relations with disease activity and joint damage. *Ann Rheum Dis*. (2001); 60: 619-26.
12. Høegh-Andersen P, Tankó LB, Andersen T, Vingsbo C, Heegaard A-M, Delaissé J-M, Christgau S. Ovariectomized Rats as a Model of Postmenopausal Osteoarthritis. Validation and Application. *Annals of Rheum Dis*. (2004); 6(2): R169-80
13. Jensen T, Hansen M, Stoltenberg M, Christgau S, Florescu A, Sommarin Y, Lorenzen I.

- Biochemical markers of connective tissue metabolism in patients with rheumatoid arthritis. Relationship to disease activity, radiographic outcome and bone mineral density. Submitted.
14. Jung M, Christgau S, Lukoschek M, Henriksen DB, Richter W. Elevated urinary concentration of collagen type II C-telopeptide fragments in patients with osteoarthritis. *Pathobiology*. (2004); 71(2): 70-6
  15. Lehmann HJ, Mouritzen U, Christgau S, Cloos PAC, Christiansen C. The effects of bisphosphonates on CartiLaps: A new marker for cartilage degradation. *Annals of Rheumatic Diseases* (2002); 61:530-533.
  16. Mazières B, Garnero P, Gueguen A, Abbal M, Berdah L, Freiburghaus C, Lequesne M, Nguyen M, Salles J-P, Vignon E, Dougados M. Molecular markers of cartilage breakdown and synovitis are strong independent predictors of structural progression of hip osteoarthritis (OA). the ECHODIAH cohort. *ACR* 2003.
  17. Mouritzen U, Christgau S, Lehmann HJ, Tankó LB, Christiansen C. CartiLaps: A novel marker of Cartilage Degradation. The influence of age, gender, menopause, hormone replacement therapy and bone mass index. *Annals Rheum Dis*. (2003); 62: 332-336.
  18. Roy-Beaudry M, Martel-Pelletier J, Pelletier J-P, M'Barek KN, Christgau S, Shipkolye F, Moldovan F. Entothelin-1 promotes osteoarthritic cartilage degradation via mmp-1 and mmp-13 induction. *Arthritis & Rheum* (2003); 48:2855-2864.

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