MutaPLATE® MTHFR (TM)

Real-Time-PCR-Kit

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14 LITERATUR

1. van der Put NM, Gabreëls F, Stevens EM, Smeitink JA, Trijbels FJ, Eskes TK, van de Heuvel LP, Blom HJ. A second common mutation in the methylenetetrahydrofolate reductase gene: an additional risk factor for neural-tube defects? *Am J Hum Genet.* 1998 May; 62(5):1044-51.
2. Khandanpour N, Willis G, Meyer FJ, Armon MP, Loke YK, Wright AJ, Finglas PM, Jennings BA. Peripheral arterial disease and methylenetetrahydrofolate reductase (MTHFR) C677T mutations: A case-control study and meta-analysis. *J Vasc Surg.* 2009 Mar; 49(3):711-8.
3. Hertfelder HJ, Gnida C, Pötzsch B, Hanfland P. MTHFR-Polymorphismus C677T: Sinn und Unsinn der Diagnostik. *Dtsch Arztebl* 2004; 101(46): A-3101 / B-2625 / C-2501. Seguí et al., *British Journal of Haematology*, 2000, 111:122-8

MutaPLATE® MTHFR (TM)

Real-Time-PCR Kit

*For the analysis of the C677T and A1298C
single nucleotide polymorphism (SNP) of the
MTHFR gene*

Valid from 2022-03-04

KF1901432
KF1901496

32/96

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1	INTENDED USE	15
2	INTRODUCTION	15
3	PRINCIPLE OF THE TEST	15
4	PACKAGE CONTENTS	15
5	EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER	16
6	TRANSPORT, STORAGE AND STABILITY	16
7	WARNINGS AND PRECAUTIONS	16
8	SAMPLE MATERIAL	17
9	REAL-TIME-PCR	17
	9.1 <i>Important points before starting</i>	17
	9.2 <i>Procedure</i>	18
	9.3 <i>Instrument settings</i>	19
10	DATA ANALYSIS	19
11	TROUBLESHOOTING	22
12	LIMITATIONS OF THE METHOD	23
13	ABBREVIATIONS AND SYMBOLS	23
14	LITERATURE	24

1 INTENDED USE

The MutaPLATE® MTHFR (TM) Real-Time PCR Kit is a molecular biological test for the parallel detection of the two common point mutations at positions 677 and 1298 of the 5,10-methylenetetrahydrofolatereductase (MTHFR) gene in open real-time PCR systems (e.g. RotorGene, SmartCycler, LightCycler, ABI, Stratagene, Amplifa). The substitution of cytosine against thymine or adenine against cytosine causes a change in the amino acid sequence from alanine to valine (Ala222Val) and in the case of A1298C from glutamine to alanine (Glu429Ala).

2 INTRODUCTION

Hyperhomocysteinemia is the result of a disturbed homocysteine metabolism often due to genetic defects. The increase of homocysteine plasma level is therefore a risk factor for cardio- or cerebrovascular complications as well as for venous thrombosis and mirgaine. Patients often carry the common C677T mutation in the MTHFR gene coding for a thermolabile variant of the MTHFR enzyme with reduced activity. The A1298C mutation in the MTHFR gene increases the homocysteine plasma level only in case of C677T heterozygosity. Due to this fact both mutations should be analysed in combination as genetically risk factors for hyperhomocysteinemia. [1], [2], [3]

3 PRINCIPLE OF THE TEST

The MutaPLATE® MTHFR (TM) Real-Time PCR Kit contains specific primers and additional material for the detection of the C677T and A1298C polymorphisms. The variable regions of the flavin adenine dinucleotide (FAD) binding site of the methylenetetrahydrofolate reductase (MTHFR) gene are amplified by PCR using genomic DNA as a template. The specific primers used in the kit flank the variable region and generate amplicons of 180 bp and 109 bp using TaqMan technology.

The standard PCR additionally contains two sequence-specific oligonucleotides labelled with a fluorescent dye (TaqMan probes). Both probes bind to the respective amplified target DNA containing the single nucleotide polymorphism (SNP). This generates a fluorescence signal which is detected by the optical unit of the real-time PCR device used. The TaqMan probe for the C and A allele (wild type) is marked with FAM (510 - 530 nm, green) and the TaqMan probe for the T and C allele (mutation) is marked with YAK (555 - 560 nm, yellow).

4 PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF1901432) or 96 (KF1901496) reactions.

Table 1: Components of the MutaPLATE® MTHFR (TM) Real-Time-PCR Kit .

Label	Lid Colour	Content	
		32	96
Enzyme mix	blue	1 x 875 µl	3 x 875 µl
Detection mix 1 C677T	yellow	1 x 368 µl	3 x 368 µl
Detection mix 2 A1298C	weiß	1 x 368 µl	3 x 368 µl
Positive control	red	1 x 30 µl	1 x 90 µl
Negative control	green	1 x 50 µl	1 x 150 µl

5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Open real-time PCR system (with plates/strips or tubes)
- Sterile PCR reaction tubes or 96-well plate/strips (white)
- Sterile reaction tubes
- Calibrated pipettes (variable volumes) and sterile disposable tips with filter
- Optional: Liquid handling system for automation

6 TRANSPORT, STORAGE AND STABILITY

The MutaPLATE® MTHFR (TM) real-time PCR kit is transported frozen on dry ice or cool packs. All components are to be stored protected from light at a minimum of -20 °C immediately after receipt. Avoid multiple freeze-thaw cycles (make aliquots if necessary). Do not use after the expiry date indicated on the package.

Be sure to protect the label detection mixes from direct sunlight during the entire test procedure.

7 WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.
- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.

- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).
- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.
- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.
- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

8 SAMPLE MATERIAL

Starting material for the MutaPLATE® MTHFR (TM) real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

9 REAL-TIME-PCR

9.1 *Important points before starting*

- Please pay attention to chapter 7 “Warnings and precautions”.
- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run all Positive Controls and one Negative Control should be included.

- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the detection mixes from exposure to light.
- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8 °C) or on ice while working.

9.2 Procedure

For amplification, two reaction tubes are needed per sample plus two additional reaction tubes per master mix for the negative and the positive controls. The following tables show the volumes to be pipetted per sample. For the analysis it is recommended to prepare a master mix for the number of samples (incl. negative and positive controls) (N) plus 10 % to compensate for inaccuracies. The master mixes are pipetted as described in tables 2 and 3:

Master mix C677T

Table 2: Preparation of master mix C677T

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 1 (yellow)	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme mix (blue)	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (green).
- For the positive control add **2 µl** of the supplied positive control (red).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding capillary.

Master mix A1298C

Table 3: Preparation of master mix A1298C

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 2 (white)	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme mix (blue)	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control (**red**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding capillary.

9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 4.

Table 4: Real-Time-PCR thermal profile

Description	Time	Temperature	Heating rate	Cycles	Acquisition
Initial Denaturation	120 s	94 °C	20 °C/s	1	none
Denaturation	10 s	94 °C	20 °C/s	45	none
Primer annealing and elongation	60 s	58 °C	20 °C/s		single
Cooling	30 s	40 °C	20 °C/s	1	-

10 DATA ANALYSIS

The results of the analysis for the **C677T** and **A1298C** polymorphisms are shown at **510-530 nm / green** and **555-560 nm / yellow** (corresponding to the Real-Time PCR instrument used).

The following figures show typical examples for both homozygous and heterozygous on a LightCycler® 2.0. With some real-time PCR devices (e.g. LightCycler® 2.0) a colour compensation file may be required.


MTHFR C667T**C-Allele** at 530 nm

Fig. 1: Evaluation of MTHFR C667T - C-Allele


MTHFR C667T**T-Allele** at 560 nm

Fig. 2: Evaluation of MTHFR C667T - T-Allele


MTHFR A1298C**A-Allele** at 530 nm

Fig. 3: Evaluation of MTHFR A1298C - A-Allele


C-Allele at 560 nm

Abb. 4: Evaluation of MTHFR A1298C - C-Allele

The Positive Control provided contains a template that is **heterozygous** for the C667T and A1298C polymorphisms.

The following three discriminations for each of the polymorphisms are possible:

C677T

1. Homozygous **C/C**:

Increase of the fluorescent signal from the **FAM** labeled TaqMan probe, no increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

2. Heterozygous **C/T**:

Increase of the fluorescent signal from the **FAM** labeled TaqMan probe and increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

3. Homozygous **T/T**:

No increase of the fluorescent signal from the **FAM** labeled TaqMan probe, increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

A1298C

1. Homozygous **A/A**:

Increase of the fluorescent signal from the **FAM** labeled TaqMan probe, no increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

2. Heterozygous **A/C**:

Increase of the fluorescent signal from the **FAM** labeled TaqMan probe and increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

3. Homozygous **C/C**:

No increase of the fluorescent signal from the **FAM** labeled TaqMan probe, increase of the fluorescent signal from the **YAK** labeled TaqMan probe.

11 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on info@immundiagnostik.com.

No fluorescence peak in the positive control or samples at about 510 - 530 nm or 555 - 560 nm.

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Detection mixes have been subjected to more than two freeze cycles or have been stored at 2-8 °C for more than four days. Repeat the analysis with a fresh aliquot or new detection mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The detection mixes were not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.

Low fluorescence peak at about 510 - 530 nm or 555-560 nm

Mix individual components carefully before use (only by pipetting several times - do not vortex!).

Cool all stock solutions appropriately during the working steps and protect the detection mixes from light irradiation.

Work on ice or with a cooled block (4 °C).

12 LIMITATIONS OF THE METHOD

The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

The accuracy of genetic tests is not 100%. However, it has been found to be over 98% accurate based on validation data for this test. Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test only analyses a selection of markers. Therefore, a negative test result of the patient does not completely exclude a risk of any kind.

13 ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleic acid		Catalog number
PCR	Polymerase chain reaction		To be used with
	Enzyme mix		Contains sufficient for <n> test
	Detection mix 1		Upper limit of temperature
	Detection mix 2		Manufacturer

Positive control

Lot number

Negative control

Consult instructions for use

Content

*In vitro* diagnostic medical device

Use by YYYY-MM-DD

14 LITERATURE

1. van der Put NM, Gabreëls F, Stevens EM, Smeitink JA, Trijbels FJ, Eskes TK, van de Heuvel LP, Blom HJ. A second common mutation in the methylenetetrahydrofolate reductase gene: an additional risk factor for neural-tube defects? *Am J Hum Genet.* 1998 May; 62(5):1044-51.
2. Khandanpour N, Willis G, Meyer FJ, Armon MP, Loke YK, Wright AJ, Finglas PM, Jennings BA. Peripheral arterial disease and methylenetetrahydrofolate reductase (MTHFR) C677T mutations: A case-control study and meta-analysis. *J Vasc Surg.* 2009 Mar; 49(3):711-8.
3. Hertfelder HJ, Gnida C, Pötzsch B, Hanfland P. MTHFR-Polymorphismus C677T: Sinn und Unsinn der Diagnostik. *Dtsch Arztebl* 2004; 101(46): A-3101 / B-2625 / C-2501. Seguí et al., *British Journal of Haematology*, 2000, 111:122-8

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com