



Data Sheet

HUMAN IgG

ANTIBODY, POLYCLONAL

Catalog no.:	G 4001
Immunogen:	Purified human IgG
Host:	Donkey
Matrix:	Serum
Specificity:	The antibody reacts specifically with human IgG
Contents:	10 ml
Known applications:	RIA (second antibody)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at:	- 20 °C
	Repeated thawing and freezing must be avoided
Last updated on:	18 August 2016

For research use only

Publishing research using G 4001? Please let us know so that we can cite your publication as a reference.

